

medicube AGE-R Booster-H Facial Treatment Device



medicube AGE-R Booster-H Facial Treatment Device User Guide

[Home](#) » [medicube](#) » medicube AGE-R Booster-H Facial Treatment Device User Guide 

Contents

- [1 medicube AGE-R Booster-H Facial Treatment Device](#)
- [2 Quick Start Guide](#)
- [3 Boosting Method](#)
- [4 Read before use!](#)
- [5 Documents / Resources](#)
 - [5.1 References](#)
- [6 Related Posts](#)

medicube⁺

medicube AGE-R Booster-H Facial Treatment Device



Quick Start Guide



STEP 01

After washing your face gently, applying the skincare product on your face.

STEP 02

Fully attached the head to the skin surface and rub the skincare product to be absorbed into the skin.

STEP 03

Refer to the boosting method, please use the device until the skincare product on your face is fully absorbed.




STEP 04

When the skincare product is fully absorbed, finish with intensive care. Bring out the natural glow of the skin with AGE-R's deep absorption technology After cleansing, apply skincare product thoroughly. Massage gently with the head of device for boosting effect in closed contact with the skin surface.



Boosting Method



Watering

Push the head of device against the skin surface and rolling it like a drawing the circle on the skin.

Boosting

Push the head of device against the skin surface and sliding it on the skin.

Focus


Push the head of device against the skin surface and pressing it one spot at a time on the skin.

Read before use!

It is recommended to start from level 1 when using it for the first time. Please adjust the level appropriately while continuously using it. 1-3 times a day


Caution

Please avoid using uvula (Adam's apple), wounds, and surgical area.



Recommended starting level

Start slowly from the lowest setting and find the appropriate level while continuously using it.




BOOSTER-H Use Cycle

You may use the Age-R Booster-H 1-3times a day. Please gradually increase the frequency of using.



Documents / Resources

	<p>medicube AGE-R Booster-H Facial Treatment Device [pdf] User Guide AGE-R Booster-H Facial Treatment Device, AGE-R Booster-H, Facial Treatment Device, Treatment Device, Device</p>
---	--

References

- [User Manual](#)